

■ 10L PLUS PILOT REACTOR, CYLINDRICAL, JACKETED, horizontal

5620

STEROGLASS

Complete reactor system on our unique open air, auto-center frame. Industrial body, high performance industrial stirring system and connections. Supplied with the components listed below and numbered on the 3D schematic drawing. Optional items are listed on the following pages.

Atex version available upon request

Horizontal version

This pilot includes:

JACKETED VESSEL DN300 10L / DISCHARGE DN50
 GLASS/PTFE VALVE DN50 WITH NO DEAD VOLUME
 COVER DN300 7 NOZZLES WITH:
 - CENTRAL DN25 FOR STIRRING,
 - DN80 PTFE NOZZLE FOR SOLIDS,
 - DN40 WITH GLASS THERMOMETER POCKET FOR PT100,
 - DN40 WITH WELL FOR PH,
 - DN40 FREE NOZZLE
 - DN25 WITH DIP PIPE AND ADAPTER FOR LIQUIDS,
 - DN50 FOR VAPOURS
 BOROSILICATE GLASS 3.3 REDUCED BEND DN100/50
 REFLUX SEPARATOR DN100
 GLASS THERMOMETER POCKET FOR PT100
 GLASS BEND 90° DN15
 5L GRADUATED FEED AND RECEIVING VESSEL
 HORIZONTAL COIL HEAT EXCHANGER 0,6MQ INC5° DN50DX
 RUPTURE DISK GRAPHITE/PFA DN50
 HIGH TORQUE-HIGH SPEED OVERHEAD STIRRER
 ASPINDLE MOTOR
 BLANTERN WITH MECHANICAL SEAL GRAPHITE/PTFE
 CSTAINLESS STEEL STIRRING SHAFT/PTFE COATED
 STAINLESS STEEL FRAME FOR REACTOR WITH COVER DN300



(r.290317)

Codice	Codice	Description	Motor type
Steroglass Fornitore			
KC-PLUS-CH10L		10L plus pilot reacelectrical	