

## ■ 20L PLUS PILOT REACTOR, CYLINDRICAL, JACKETED, horizontal

5623

### STEROGLASS

Complete reactor system on our unique open air, auto-center frame. Industrial body, high performance industrial stirring system and connections. Supplied with the components listed below and numbered on the 3D schematic drawing. Optional items are listed on the following pages.

horizontal version  
Atex version available upon request

This pilot includes:

JACKETED VESSEL DN300 20L / DISCHARGE DN50  
GLASS/PTFE VALVE DN50 WITH NO DEAD VOLUME  
COVER DN300 7 NOZZLES WITH:  
- CENTRAL DN25 FOR STIRRING,  
- DN80 PTFE NOZZLE FOR SOLIDS,  
- DN40 WITH GLASS THERMOMETER POCKET FOR PT100,  
- DN40 WITH WELL FOR PH,  
- DN40 FREE NOZZLE  
- DN25 WITH DIP PIPE AND ADAPTER FOR LIQUIDS,  
- DN50 FOR VAPOURS  
BOROSILICATE GLASS 3.3 REDUCED BEND DN100/50  
REFLUX SEPARATOR DN100  
GLASS THERMOMETER POCKET FOR PT100  
GLASS BEND 90° DN15  
n° 2 5L GRADUATED FEED AND RECEIVING VESSEL  
HORIZONTAL COIL HEAT EXCHANGER 0,6MQ INC5° DN50DX  
RUPTURE DISK GRAPHITE/PFA DN50  
HIGH TORQUE-HIGH SPEED OVERHEAD STIRRER  
A SPINDLE MOTOR  
B LANTERN WITH MECHANICAL SEAL GRAPHITE/PTFE  
C STAINLESS STEEL STIRRING SHAFT/PTFE COATED  
STAINLESS STEEL FRAME FOR REACTOR WITH COVER DN300

(r.290317)

Codice	Codice	Description	Motor type
Steroglass Fornitore			
KC-PLUS-CH20L		20L plus pilot	Electrical

