

■ 50L PLUS PILOT REACTOR, CYLINDRICAL, JACKETED DN450 horizontal

5630

STEROGLASS

Complete reactor system on our unique open air, auto-center frame. Industrial body, high performance industrial stirring system and connections. Supplied with the components listed below

-Atex version available upon request
horizontal version

This pilot includes:

50L DISMOUNTABLE JACKETED VESSEL DN450/400 50L DISCHARGE DN50
GLASS/PTFE VALVE DN50 WITH NO DEAD VOLUME
COVER DN300 7 NOZZLES WITH:
- CENTRAL DN25 FOR STIRRING,
- DN80 PTFE NOZZLE FOR SOLIDS,
- DN40 WITH GLASS THERMOMETER POCKET FOR PT100,
- DN40 WITH WELL FOR PH,
- DN40 FREE NOZZLE
- DN25 WITH DIP PIPE AND ADAPTER FOR LIQUIDS,
- DN50 FOR VAPOURS
BOROSILICATE GLASS 3.3 REDUCED BEND DN100/50
REFLUX SEPARATOR DN100
GLASS THERMOMETER POCKET FOR PT100
GLASS BEND 90° DN15
5L GRADUATED FEED AND RECEIVING VESSEL
HORIZONTAL COIL HEAT EXCHANGER 1,5MQ DN80
RUPTURE DISK GRAPHITE/PFA DN50
HIGH TORQUE-HIGH SPEED OVERHEAD STIRRER
ASPINDLE MOTOR
BLANTERN WITH MECHANICAL SEAL GRAPHITE/PTFE
CSTAINLESS STEEL STIRRING SHAFT/PTFE COATED
STAINLESS STEEL FRAME FOR REACTOR WITH COVER DN300

(r.300317)



Codice	Codice	Description	Motor type
Steroglass Fornitore			
KC-PLUS-		50L plus pilot	Electrical