

Name	Value
Pharmacy refrigerator	MD40DC-P600BEQMXXXXX1SN4B
EAN code	4770104040061
Gross capacity, l	385
Net capacity, l	290
Factory temperature setting, °C	+5
Fluctuation, °C	0,5
Gradient, °C	2,9
Increase in temperature from +5 °C to +10 °C (empty, +25 °C TU)	40 min.
Decrease in temperature from +25 °C to +5 °C (empty, +25 °C TU)	40 min.
Recovery time after 1 min of door opening (empty, +25 °C TU, according to EN 60068-3)	6 min.
Positive temperature deviation from the setpoint	1.3
Negative temperature deviation from the setpoint	-1.2
Heat output, kJ/h	252
Ambient temperature range, °C	from +10 to +35
Airborne acoustical noise emissions, dB(A) re1 pW	44
Refrigerant	R600a
Daily consumption, kWh/24h	0.953
Annual energy consumption, kWh/year	348
Rated voltage, V	220 - 240
Rated frequency, Hz	50
Self closing door	Yes
Access port	Yes
Cooling system	Fan forced cooling
Defrost	Automatic
Volt-free alarm contact	No
Interior light	Vertical LED light strip on right side
Malfunction warning signal	Optical and acoustical
Protection against low temperature	Yes
Number of shelves	5 adjustable +1 bottom
Shelf size (W x D), mm	380 x 495, bottom 445 x 180
Shelf installation step, mm	12.5
Maximum load per shelf, kg	40
Shelf material	Wire shelves coated by white polymer paint
Side wall material	Steel, white painted
Inner cabinet material	HIPS polystyrene
Type of door	Full-panel door

Name	Value
Type of lock	Mechanical
Adjustable foot at the front	2
Wheels behind	2
External dimensions (H x W x D), mm	2025 x 600 x 620
Internal dimensions (H x W x D), mm	1645 x 500 x 418
Package dimensions (H x W x D), mm	2080 x 620 x 630
Weight net, kg	76.5
Weight gross, kg	82
Pharmacy refrigerator packed on a pallet	MD40DC-P600BEQMXXXXX1SN4A
EAN code	4770104040078
Packaged on pallet dimensions (H x W x D), mm	2210 x 620 x 710
Weight gross (packaged on pallet), kg	90,5